



CATALOGUE

smartbone® **smart**screws™ **smart**tools™ **smart**kit™

smartbone® Block

ITEM	SIZE	Q.TY
SMB011005	7 X 7 X 7 mm	1
SMB011010	10 X 10 X 10 mm	1
SMB011020	10 x 10 x 20 mm	1
SMB011030	10 x 20 x 20 mm	1
SMB011110	14 x 12 x 6 mm	1
SMB011130	14 x 12 x 8 mm	1
SMB011160	14 x 12 x 12 mm	1
SMB011180	14 x 12 x 20 mm	1
SMB011190	14 x 12 x 24 mm	1
SMB011210	15 x 30 x 20 mm	1
SMB011220	15 x 30 x 30 mm	1
SMB011230	15 x 30 x 40 mm	1
SMB011250	15 x 30 x 60 mm	1
SMB011310	16 x 14 x 6 mm	1
SMB011330	16 x 14 x 8 mm	1



smartbone® Microchips

ITEM	SIZE	Q.TY
SMC010050	0.25 - 1 mm	0.5 cc (~ 0.3 g)
SMC010100	0.25 - 1 mm	1 cc (~ 0.6 g)
SMC010110	0.25 - 1 mm	2 cc (~ 1.2 g)
SMC010190	0.25 - 1 mm	5 cc (~ 3 g)
SMC010221	1 - 2 mm	1 cc (~ 0.38 g)
SMC010200	1 - 2 mm	2 cc (~ 0.75 g)
SMC010201	1 - 2 mm	4 cc (~ 1.5 g)



smartbone® Granules

ITEM	SIZE	Q.TY
SMC010300	2 - 4 mm	5 cc (~ 1.6 g)



smartbone® Plate

ITEM	SIZE	Q.TY
SMP013040	4 x 10 x 10 mm	1
SMC013010	3 x 25 x 15 mm	1



smartbone® Special U Shape

ITEM	SIZE	Q.TY
SMS016010	15 X 6 X 4 mm	1



smartbone® Special L Shape

ITEM	SIZE	Q.TY
SMS016210	15 X 7 X 6 mm	1



smartbone® Special C Shape

ITEM	SIZE	Q.TY
SMS016100	ø 8 x h 12 mm	1



smartbone® Cylinder

ITEM	SIZE	Q.TY
SMC012010	ø 11 x h 20 mm	1
SMC012030	ø 14 x h 20 mm	1
SMC012050	ø 11 x h 10 mm	1
SMC012060	ø 14 x h 14 mm	1



smartbone® Rod

ITEM	SIZE	Q.TY
SMR015010	35 x 4 x 3 mm	1





smartscrews™

smartscrews™ - Titanium

ITEM	DESCRIPTION	SIZE	Q.TY
SMK0715T	1 package of 5 SmartScrews™	Lenght 7 mm / ø 1,5 mm TI ASTM Gr5	1
SMK1115T	1 package of 5 SmartScrews™	Lenght 11 mm / ø 1,5 mm TI ASTM Gr5	1
SMK1515T	1 package of 5 SmartScrews™	Lenght 15 mm / ø 1,5 mm TI ASTM Gr5	1

smartscrews™ - Steel

ITEM	DESCRIPTION	SIZE	Q.TY
SMK075S	1 package of 5 SmartScrews™	Lenght 7 mm / ø 1,5 mm AISI 316L	1
SMK1115S	1 package of 5 SmartScrews™	Lenght 11 mm / ø 1,5 mm AISI 316L	1
SMK1515S	1 package of 5 SmartScrews™	Lenght 15 mm / ø 1,5 mm AISI 316L	1






smarttools™

Drill

ITEM	DESCRIPTION	Q.TY
SMD08	osteosynthesis drill ø 0,8 mm - h 6 mm (AISI630)	1
SMD12	osteosynthesis drill ø 1,2 mm - h 12 mm (AISI630)	1
SMD15	osteosynthesis drill ø 1,5 mm - h 12 mm (AISI630)	1
SMD25	round burr ø 2,5 mm (M340)	1
SMD23	shaping bone drill ø 2,3 mm (M340)	1



Screwdriver

ITEM	DESCRIPTION	Q.TY
SCAR	implant screwdriver 22 mm	1
SCCA	contra-angle connector Lenght 25,3 mm	1
SCDC	screwdriver for handpiece Lenght 25,3 mm	1
SCDM	manual screwdriver Lenght 160 mm	1





smarttools™

Bone Cutter

ITEM	SIZE	Q.TY
SCTR	Lenght 155 mm (mod. Ruskin mini)	1



Surgical Probe

ITEM	SIZE	Q.TY
SCSN	Lenght 55 mm (mod. University of North Carolina)	1



smartkit™

Kit smarttools™

ITEM	DESCRIPTION	Q.TY
SCKT	aluminum tray with sliding lid	1



Osteosynthesis smarttools™

ITEM	DESCRIPTION	Q.TY
IBI-KIT01	<ul style="list-style-type: none"> - alluminum tray with sliding lid - manual screwdriver - screwdriver for handpiece - implant screwdriver - contra-angle connector - osteosynthesis drill \varnothing 0,8 mm - h 6 mm - osteosynthesis drill \varnothing 1,2 mm - h 12 mm - osteosynthesis drill \varnothing 1,5 mm - h 12 mm - round burr \varnothing 2,5 mm - shaping bone drill \varnothing 2,3 mm 	1





Via Cantonale 67
CH-6805 Mezzovico-Vira
Switzerland

t. +41 91 93 06 640
f. +41 91 220 72 00
www.ibi-sa.com
info@ibi-sa.com

VAT code: CHE-112-388-119 IVA